

# OPTICLUX 10



medicinelaser<sup>®</sup>  
LATAM



i NOVO<sup>®</sup>  
AWARDS  
2017



reddot winner 2018

## Details shown in the best light. OPTICLUX 10

The OPTICLUX 10 magnifier light combines a robust magnifier with outstanding light quality and ergonomic operation. As a certified medical product, it is used in dermatology, in clinics and hospitals for human, dental and veterinary medicine, in laboratories and veterinary medicine, in laboratories, but also for central steril supply department (CSSD).



OPTICLUX 10	technical data
<b>luminous power *</b>	6,000 Lux @ 0.15 m
<b>colour temperature *</b>	<ul style="list-style-type: none"> <li>• 6500 K</li> <li>• optional: 365 nm (Wood's Light)</li> </ul>
<b>lens dimension</b>	Ø 16 cm
<b>colour rendering index</b>	Ra > 95
<b>usage</b>	membrane switch on light head
<b>class of protection</b>	<ul style="list-style-type: none"> <li>• version with cover: I</li> <li>• version with pin: II</li> </ul>
<b>protection rating</b>	IP 20

special features	
<b>OPTICLUX 10</b>	<ul style="list-style-type: none"> <li>• dimmable</li> <li>• maintenance-free</li> <li>• Wood's Light mode (optional)</li> <li>• lens cover</li> <li>• 1.9x magnification</li> </ul>
<b>variations</b>	Further connector variants are available on request.
<b>accessories</b>	<ul style="list-style-type: none"> <li>• table clamp</li> <li>• universal mount</li> <li>• fasting hardware</li> <li>• wall mount</li> <li>• roller stand</li> </ul>

# OPTICLUX 10 – variations



OPTICLUX 10 C  
ceiling version



OPTICLUX 10 P  
pin version



## Product highlights



**Optimally protected**  
sturdy plastic cover  
protects the lens from  
scratches and sunlight




**360° turnover axis**  
for maximum light  
availability




**Especially for dermatology**  
two illumination modes in  
one device:  
optional Wood's Light for  
examination of  
fluorescent fungal or  
bacterial infection on the  
skin or scalp and  
pigmentary changes.

**Detailed description**  
to Wood's Light on page 110.



**Ergonomic handling**  
spring arm and internal  
friction joints allow  
comfortable handling  
without any risk of the  
light head drifting



**Flexible positioning**  
long lamp arm makes a  
large radius of action up  
to 142 cm possible

**medicinelaser**  
LATAM  
[www.medicinelaser.com](http://www.medicinelaser.com)



ARGENTINA  
Av. Doctor Honorio  
Pueyrredón 707, Caballito,  
CABA, Argentina.  
M: +54 11 3404 7158

CHILE  
Sucre 1355, Ñuñoa,  
Santiago, Chile.  
O: +56 2 2269 1981  
M: +56 9 4409 7105

ECUADOR  
Av. República 1058,  
Edificio La Fontana,  
Quito, Ecuador.  
M: +593 99 120 4597

PERÚ  
Calle Almirante Lord  
Cochrane 322,  
Miraflores, Lima, Perú.  
O: + 51 1 704 80 14  
M: +51 983 552 813

[ventas@medicinelaser.com](mailto:ventas@medicinelaser.com)