

TRIANGO 80

medicinelaser[®]
LATAM



reddot winner 2012

A proven classic. TRIANGO 80

With 80,000 lux at a distance of one meter and a light field diameter of 18 cm, the medical staff has a powerful treatment light at its disposal.

The illumination field is pleasantly homogeneous and continuously dimmable.

The arm system enables simple and accurate alignment of the lamp head position.

The handle makes it easy for right-handed and left-handed users to position the light.



TRIANGO 80	technical data
luminous power *	80,000 Lux @ 1.0 m
colour temperature *	ca. 4500 K
lighting field @ 1.0 m *	Ø 18 cm
depth of illumination	124 cm
colour rendering index	Ra > 95
usage	rotary switch on the luminaire head
class of protection	I
protection rating	IP 43 (luminaire head) IP 20 (UPSmed pack)

special features

TRIANGO 80-1	<ul style="list-style-type: none"> continuously dimmable
fastening	<ul style="list-style-type: none"> ceiling wall roller stand roller stand UPSmed pack

TRIANGO 80 – variations



TRIANGO 80 C ceiling version



TRIANGO 80 W wall mount



D

TRIANGO 80 F roller stand



TRIANGO 80 F Battery-UPSmed Pack

