

VISIANO 10



medicinelaser*
LATAM 



Reduced to the essentials. VISIANO 10

The VISIANO 10 is an examination light that is reduced to the essentials, but still delivers professional lighting technology. The illuminance, with 50,000 lux at a distance of 0.5 m, is suitable for general examinations but also for intensive care and neonatology. The ergonomic handle makes it easy to align the light field.



VISIANO 10	technical data
luminous power *	50,000 Lux @ 0.5 m
colour temperature *	4400 K
lighting field @ 1.0 m *	Ø 18 cm
colour rendering index	Ra > 93
usage	rocker switch (I/O) on arm
class of protection	II
protection rating	IP 20


Special features	
VISIANO 10-1 P S10	<ul style="list-style-type: none"> • handle • luminaire head can be swivelled 280°
fastening	<ul style="list-style-type: none"> • pin universal
accessories	<ul style="list-style-type: none"> • wall mount • roller stand • fastening hardware • universal mount • table clamp

VISIANO 10 – variations



VISIANO 10 P

pin version | with accessory roller stand



Various connections

The VISIANO 10 can be combined with Derungs accessories. Depending on the requirements, it can be fixed directly to the wall with a wall mount or a rail mount. The Derungs roller stand is suitable for mobile use.

Easy cleaning

Thanks to its closed and minimalistic design, the VISIANO 10 can be cleaned in a time-saving manner. Dust and liquids cannot enter the interior. This means that nothing stands in the way of long-lasting use.

Product highlights



Economical
due to low energy
consumption and
high-quality components

Durable flex arm
for safe positioning and
durable use

Large radius of action
thanks to 280° rotatable
lamp head and special
arm

medicinelaser
LATAM
www.medicinelaser.com



ARGENTINA
Av. Doctor Honorio
Pueyrredón 707, Caballito,
CABA, Argentina.
M: +54 11 3404 7158

CHILE
Sucre 1355, Ñuñoa,
Santiago, Chile.
O: +56 2 2269 1981
M: +56 9 4409 7105

ECUADOR
Av. República 1058,
Edificio La Fontana,
Quito, Ecuador.
M: +593 99 120 4597

PERÚ
Calle Almirante Lord
Cochrane 322,
Miraflores, Lima, Perú.
O: + 51 1 704 80 14
M: +51 983 552 813

+ OFICINAS CENTRALES

ventas@medicinelaser.com