

VISIANO 20



medicinelaser^{*}
LATAM



reddot winner 2014

Outstanding in design and function. VISIANO 20

The VISIANO 20 is the premium product for general examinations.

The light output of 60,000 lux at a distance of 0.5 m is outstanding in this class.

This makes it suitable for versatile applications: general medicine, dermatology, gynecology, urology, neonatology and anesthesia. But also in intensive care, emergency and in the preparation for surgery, it provides excellent illumination and LEDs in two different light colours.



VISIANO 20	technical data
luminous power *	60,000 Lux @ 0.5 m
colour temperature *	3500 K & 4500 K
lighting field @ 1.0 m *	Ø 21 cm
colour rendering index	Ra > 95
usage	touch panel on light head
class of protection	ceiling version: I pin version: II
protection rating	IP 20

special features	
VISIANO 20-2	<ul style="list-style-type: none"> • 2 colour temperatures • dimmable
fastening	<ul style="list-style-type: none"> • ceiling • pin universal
accessories	<ul style="list-style-type: none"> • wall mount • roller stand • fastening hardware • universal mount • table clamp

VISIANO 20 – variations



VISIANO 20 C ceiling version



VISIANO 20 P pin version



Versatile use

The mounting options leave nothing to be desired:
The VISIANO Pin version can be combined with versatile accessories e.g. roller stand, rail mount, wall mount.
The ceiling version rounds off the mounting options.

High quality lighting technology

VISIANO 20 is equipped with LEDs in two colour temperatures as standard. The special optics provide a glare-free and homogeneous illumination field.
High-quality materials and long-life LEDs ensure a flawless use.

Product highlights



Colour change
for optimal contrast vision



Unique diamond optic
for best colour homogeneity of the light



360° turnover axis
for maximum light availability



Spring-relieved friction joint
for secure hold

medicinelaser
LATAM

www.medicinelaser.com



ARGENTINA
Av. Doctor Honorio
Pueyrredón 707, Caballito,
CABA, Argentina.
M: +54 11 3404 7158

CHILE
Sucre 1355, Ñuñoa,
Santiago, Chile.
O: +56 2 2269 1981
M: +56 9 4409 7105

ECUADOR
Av. República 1058,
Edificio La Fontana,
Quito, Ecuador.
M: +593 99 120 4597

PERÚ
Calle Almirante Lord
Cochrane 322,
Miraflores, Lima, Perú.
O: + 51 1 704 80 14
M: +51 983 552 813

+ OFICINAS CENTRALES

ventas@medicinelaser.com